

# Pro Water Baths

RS-SP-100 & RS-PB-50



## Revolutionary Science

*Manufacturer of Precision Laboratory Equipment*

Quick, durable, compact! The modest sized Mini Pro Bath is perfect for the individual or applications that don't require a large bath. The high impact polypropylene vented cover allows for tall beakers, bottles or flasks to saturate in a precision temperature controlled environment. Two year warranty included.

### Mini ProBath RS-PB-50 specifications

Precision:	Internal temperature accuracy of $\pm 2.0^{\circ}\text{C}$
Temperature Range:	Ambient to $99.9^{\circ}\text{C}$
Capacity:	0.5L
Net Weight:	1kg (2.3lbs)
Power Requirements:	120 VAC, 50-60 Hz, 200 W, grounded cord
Durability:	Rugged exterior is resistant to dents and cracks
Portability:	Great for transportation
Dimensions:	
Exterior:	20.3cm diameter x 11.4cm tall 8" diameter x 4.5" tall
Interior:	10.5cm diameter x 6.5cm tall 4.1" diameter x 2.6" tall



With a powerful, precision controlled 1250 watt heating element, the Speed Pro Bath is one of the fastest water baths on the market today. Built from medical grade, hand polished stainless steel, it not only looks sharp, but can also be used for countless applications in the laboratory setting. Includes an unbreakable polycarbonate cover. Two year warranty included.

### Speed Pro RS-SP-100 specifications

Precision:	Internal temperature uniformity of $\pm 0.4^{\circ}\text{C}$
Display:	3 digit, $0.1^{\circ}\text{C}$ increments
Speed:	Temperature reaches boiling in about 10 minutes
Temperature Range:	Ambient to $99.9^{\circ}\text{C}$
Capacity:	7.9L
Net Weight:	4kg (9lbs)
Power Requirements:	100-240 VAC, 50-60 Hz, 1250 W, grounded cord
Durability:	Stainless steel chamber, rugged ABS exterior
Dimensions:	
Exterior:	38cm wide x 30.5cm deep x 25.4cm tall 15" wide x 12" deep x 10" tall
Interior:	23cm diameter x 19cm tall



## Revolutionary Science

17319 Lake Blvd.  
Shafer, MN 55074  
Sales: (651) 674-3109  
Fax: (775) 218-7015

E-mail: [Sales@revsci.com](mailto:Sales@revsci.com)  
Website: [www.revsci.com](http://www.revsci.com)

**ISO 9001**  
REGISTERED FIRM