

# Electrophoresis RS-PS-75

## Power Supply



### Revolutionary Science

*Manufacturer of Precision Laboratory Equipment*

The Revolutionary Science power supply is the ideal power supply for running routine electrophoresis gels. The RS-PS-75 is capable of running two gels simultaneously at the same or different voltages. It is equipped with a red warning LED light that illuminates when output current exceed 350mA. Rapid discharge feature eliminates the possibility of electric shock when the unit is turned off. Two year warranty included.

### RS-PS-75 specifications

Max Power:	52.5 watts
Max Current:	350mA
Output Voltages:	0 volts, 15 volts, 75 volts, 150 volts
Output Jack Sets:	2
Fuse:	Glass encapsulated
Net Weight:	0.68kg (1.5lbs)
Simplicity:	Easy to use slide switches
Dimensions:	7cm tall x 12.7cm wide x 10.2cm deep 2.8" tall x 5" wide x 4" deep



### Revolutionary Science

17319 Lake Blvd.  
Shafer, MN 55074  
Sales: (651) 674-3109  
Fax: (775) 218-7015

E-mail: [Sales@revsci.com](mailto:Sales@revsci.com)  
Website: [www.revsci.com](http://www.revsci.com)

**ISO 9001**  
REGISTERED FIRM